

# Malaria

The following tests are available:

- **Malaria antigen, thick drop, microscopy/immunoassay<sup>^</sup>**

Indication: Suspicion of malaria, acquired fever in areas with malaria risk

Material: EDTA blood

Preanalytics: collection during fever peaks

Stability: 3 days at 2 to 8°C

TAT: same day, FML

Method: Blood smear, thick drop, IA

- **Malaria IgG antibodies**

Indication: Clarification of earlier fever attacks after traveling in malaria risk areas, test is not suitable for the diagnosis of acute malaria.

Material: 1 ml serum

TAT: 7-10 days \*

Method: EIA

Units: Index

Ref.- range: see report

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit

<http://www.fml-dubai.com/parameter-listings/>