

# Salmonella

The following tests are available:

- **Salmonella antibodies (Widal reaction)**

Indication: screening test for post-infectious salmonella disorder, clarification of arthralgia, aneurysms

Material: 1 ml serum

TAT: same day, FML

Method: Widal reaction

Units: qualitative

Ref.- range: negative

- **Salmonellae in stool, bacteriological determination**

Indication: diarrheas (e.g. aqueous or bloody)

Material: 5 g stool

TAT: 2 - 3 days , FML

Method: cultivation on selective media, serum typing of O and H phases by slide agglutination

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit <http://www.fml-dubai.com/parameter-listings/>