

TSH thyroid stimulating hormone

Material: 1 ml serum

Stability: 14 days at 2 to 8°C

TAT: same day, FML

Method: ECLIA

Unit: mIU/L

Ref.- range: 0.30 - 4.20

Note: If the patient is taking multivitamins or dietary supplements containing high dose of Biotin (> 5 mg), the patient should stop taking it for at least 24 hours , before having the blood collection.

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit
<http://www.fml-dubai.com/parameter-listings/>