

# Alcohol

The following tests are available:

- **Alcohol in serum**

Material: 2ml serum

Preanalytics: do not use alcohol pads for skin disinfection, the collection tube must be closed until analysis (to avoid evaporation).

Stability: 14 days at 2-8°C

TAT: same day, FML

Method: photometry

Units: g/l

Ref.-range: <0.1

Note: Not recommended for forensic purposes, only for medical indications.

- **Alcohol in urine**

Material: 10 ml urine

Stability: 30 days at 2-8°C

TAT: same day, FML

Method: photometry

Units: g/l

Ref.-range: <0.1

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit <http://www.fml-dubai.com/parameter-listings/>