

Crosslinks

see also **Procollagen**

General:

Pyridinolines are resorption products of collagen and reflect the quantity of degraded collagen. They are considered markers for bone (deoxypyridinoline) and cartilage (pyridinoline) resorption.

Indication: Increased reorganization / resorption / reduction of bone

Material: 10 ml urine

Preanalytics: 2nd morning urine

TAT: 7-10 days*

Method: HPLC

Units: µg/l

Ref.- range: see report

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit <http://www.fml-dubai.com/parameter-listings/>