

Ethanol Ethyl Alcohol C₂H₅OH

The following tests are available:

- **Alcohol in serum**

Preanalytics: do not use alcohol patches for skin disinfection, the collection tube must be closed until investigation (to avoid evaporation).

Stability: 14 days at 2 to 8°C

Material: 1 ml serum

TAT: same day, FML

Method: photometry

Units: g/l

Ref.- range: <0.1

Note: Not recommended for forensic purposes, only for medical indications.

- **Alcohol in urine**

Material: 1 ml urine

TAT: same day, FML

Method: photometry

Units: g/l

Ref.- range: <0.1

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit
<http://www.fml-dubai.com/parameter-listings/>