

Triglycerides

- **Triglycerides in serum**

Material: 1 ml serum

Stability: 10 days at 2 to 8°C

TAT: same day, FML

Preanalytics: blood collection after at least 12-hour period of fasting.

Method: photometry

Ref.- range: see report

- **Triglycerides in aspirate[^]**

Material: 1 ml aspirate

TAT: same day, FML

Method: photometry

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit <http://www.fml-dubai.com/parameter-listings/>