

Freiburg Medical Laboratory ME LLC, P.O.Box 3068, Dubai

Tel: 04 396 2227

Fax: 04 396 2228



# Anti mitochondrial antibodies AMA

#### General:

Anti-mitochondrial antibodies (AMA) are antibodies formed against mitochondria, primarily in cells of the liver. AMAs actually bind to protein antigens within the inner lining of the mitochondria. The presence of AMAs in the blood or serum is indicative of several autoimmune diseases such as primary biliary cirrhosis. It is present in about 95% of cases. Primary biliary cirrhosis is seen primarily in middle-aged women, and in those afflicted with other autoimmune diseases. PBC is an autoimmune disorder, a condition in which the immune system attacks the body's own cells, or in this case parts of the cells. It is postulated that xenobi-otically-induced and/or oxidative modification of mitochondrial autoantigens is a critical step leading to loss of tolerance causing AMAs. In acute liver failure AMAs are found against all major liver antigens. A subtyping of the antimitochondrial antibodies (AMA) is recommended in cases of positive screening. The most important subtypes are AMA-M2, AMA-M4, AMA-M9.

The following tests are available:

## AMA Screening

Material: 1 ml serum

Stability: up to 14 days at 2 to 8°C

TAT: 3 days, FML

Method: IFT

## AMA-M2

Material: 1 ml serum

Stability: up to 14 days at 2 to 8°C

TAT: 3 days, FML

Method: dot blot

#### AMA-M4

Material: 1 ml serum TAT: 7-10 days\*

Method: blot

### AMA-M9

Material: 1 ml serum

Page 1 of 2 Updated 27/03/2022





Copyright Freiburg Medical Laboratory

Freiburg Medical Laboratory ME LLC, P.O.Box 3068, Dubai

Tel: 04 396 2227

Fax: 04 396 2228

TAT: 7-10 days\*

Method: blot

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit <a href="http://www.fml-dubai.com/parameter-listings/">http://www.fml-dubai.com/parameter-listings/</a>

Page 2 of 2 Updated 27/03/2022

