

ص.ب: ٢٠٦٨، دبي - الإمارات العربية المتحدة هاتف : ۲۲۲۷ ۲۹۲ ۶ و فاکس: ۲۲۲۸ ۲۹۲ ۶۰

P.O. Box. 3068, Dubai - U.A.E. Tel. 04 396 2227 Fax. 04 396 2228 E-mail, info@fml-dubai.com Website www.fml-dubai.com

مخنبر فرايبورغ الطبي الشرق الأوسط (د.م.م) FREIBURG MEDICAL LABORATORY MIDDLE EAST (L.L.C.)

رأس المال ٦٠٠،٠٠٠ درهم مدقوع بالكامل – سجل تجاري ٩٩٥٥ (١٩٤٥ - ١٩٥٥ Capital: AED 600,000 Fully Paid

A member of SYNLAB

Semen Analysis Report

Patient Information		Laboratory Information				
Name	Sample Report	Ordering Physician	PD Dr. med. habil. M. Jaksch			
Your ID		Clinic Name	Freiburg Medical Laboratory			
Patient ID		14				
Birth Date	01.01.1981	Release Date	13.11.2018			
Abstinence (# of days)	4					

Sample Processing			
Date Sample Collected	13.11.2018	Time Sample Collected	14:25
Date Sample Received	13.11.2018	Time Sample Received	14:33
Date Test Performed	13.11.2018		

Sample Collection and Initial Examination			Test Results		Ref. Value			
		SQA-V Automated Test			Laboratory	tory WHO 5th	↑ ↓	
Sample/Accession #	A000000	Volume	1.5	ml		>/= 1.5		
Color	Grevish white	Viscosity	Normal					
Liquefaction	1 hour	рН	8.0			>=7.2		
Appearance	Homogeneous	WBC Conc.	<1	M/ml		<1	- 7	

	ж.				Te	st Re	esults						
			Ref. Value		T				Ref. Value				
SQA-V Automated Test		Test Results		Manufacturer	WHO 5th >= 15.0	↑↓	Manual Test Agglutination	Test Results		Laboratory	WHO 5th	↑↓	
Sperm Conc.		101.3 M/m						-					
Total Motility (PR+NP)		42	%		>= 40		Aggregation	-					
Prog. Motility (PR)		34	%		>= 32		Round cells	N.A.	M/ml				
Nonprog. Motility (NP)		8	%				RBC	N.A.	M/ml				
Immotility (IM)		58	%				Vitality	51	%		>/= 55	. ↓	
Morph. Normal Forms*	WHO 5th	8	%		>= 4		Manual Morphology						
MSC		42.0	M/ml	>=6			% Normal Forms	55	%	>30			
PMSC		34.7	M/ml	>=5			Head Defects	35	%	- 30		-	
FSC		6.5	M/ml				Neck/Midpiece	8	%				
Velocity		42	mic/sec	>=5			Tail Defects	2	%				
SMI		202					Cytoplasmic Droplets	0	%				
Sperm # (TSC)		151.9	M		>= 39.0		Pinheads	0	%				
Motile Sperm		63.1	М	>=16									
Prog. Motile Sperm		52.1	М	>=12			WBC	3-5	/hpf	0-5			
Functional Sperm		9.8	М				Bacteria	+					
Morph. Normal Sperm		12.2	М	>=2							-		

Interpretation

Normal sperm concentration and acceptable range of morphological normal forms.

Comment

Please note that the newly introduced SQA-V automated analyzer results are not necessarily compatible with the manual test results. Please also note that all results should be related to the WHO (5th) edition.

TSC: Total Sperm Concentration (Count)

PMSC: Progressive Motile Sperm Concentration

SMI: Sperm Motility Index

MSC: Motile Sperm Concentration

FSC: Functional Sperm Concentration

Technical Validation by Medical Technologist (Supervisor of the Department)

Dr. Nehmat ElBanna Specialist Clinical Pathology (U/S) Licence No: DHA-P-0084548



PD Dr. med. habil. M. Jaksch **Associate Professor Medical Director** Licence No: DHA-LS-240710 Page 1/1