

Blood differential manual

Material: EDTA-blood, fresh or smear.

Stability: Max 5 hours at 20 to 25°C

TAT: Same day, FML

Ref.- range: Neutrophils 40 - 70 % (Adults)

Band forms 0 - 5 %

Eosinophils 0 - 4 %

Basophils 0 - 2 %

Lymphocytes 20 - 40 %

Monocytes 4 - 12 %

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit
<http://www.fml-dubai.com/parameter-listings/>